



TEST: Anti- β_2 -Glycoprotein I IgG, IgM or IgA Antibody Detection in Human Serum by ELISA

PRINCIPLE:

β_2 glycoprotein I (β_2 GPI) is a 50-kDa glycoprotein that is present in the plasma and binds as a cofactor to anionic phospholipids (such as cardiolipin). It also binds to the cell membrane of activated platelets and endothelial cells which are involved in the antiphospholipid syndrome. During pregnancy, β_2 GPI can bind to trophoblast cells. The subsequent binding of anti- β_2 GPI antibodies may cause the decrease of hCG secretion from trophoblast cells which can lead to defective placentation. β_2 GPI also functions as a natural anti-coagulant. Therefore, patients with anti- β_2 GPI antibodies may be at increased risk for thrombosis. The anti- β_2 GPI assay may be used in conjunction with the anti-phospholipid antibodies assay for aiding in the diagnosis of antiphospholipid syndrome. The three isotypes of anti- β_2 GPI (IgG, IgM and IgA) are measured in human serum by ELISA.

Purified β_2 GPI antigen is bound to the wells of a polystyrene microwell plate under conditions that will preserve the antigen in its native state. Pre-diluted controls and diluted patient sera are added to separate wells, allowing any β_2 GPI IgG antibodies present to bind to the immobilized antigen. Unbound sample is washed away and an enzyme labeled anti-human IgG conjugate is added to each well. A second incubation allows the enzyme labeled anti-human IgG to bind to any patient antibodies, which have become attached to the microwells. After washing away any unbound enzyme labeled anti-human IgG, the remaining enzyme activity is measured by adding a chromogenic substrate and measuring the intensity of the color that develops. The assay can be evaluated spectrophotometrically by measuring and comparing the color intensity that develops in the patient wells with that of a five point calibration curve.

SPECIMEN REQUIREMENTS:

2ml collected in a serum separator tube (gel barrier). Separate serum from cells ASAP or within 2 hours of collection by centrifugation. Stability after separation from cells: Ambient: 8 Hours; Refrigerated: 48 Hours; Frozen: 1 year (avoid repeated freeze/thaw cycles).

REJECTION CRITERIA:

Plasma or other body fluids. Gross hemolysis

METHOD: Enzyme Linked Immunoassay (ELISA)

REFERENCES:

1. H. Bas de Laat et al. *Clin Immunol* 2004;112:161-168
2. Di Simone et al. *Ann Rheum Dis* 2005;64:462-467

Normal Range: Negative.

Results are reported as negative or positive.

Turnaround Time: 7-10 business days