

TEST: VARICELLA-ZOSTER IgG

PRINCIPLE:

Varicella (chickenpox) and zoster (shingles) represent different clinical manifestations of infection with the same agent, Varicella-zoster virus (VZV), a member of the *Herpesviridae*.

Varicella occurs most frequently in children and is characterized by a generalized vesicular exanthema often accompanied by fever. Zoster usually occurs in adults or immunocompromised patients (including those with AIDS) and consists of painful, circumscribed eruption of vesicular lesions with accompanying inflammation of associated dorsal root or cranial nerve sensory ganglia. Varicella is the primary infection with VZV, whereas zoster is a secondary infection due to reactivation of latent VZV sensory ganglia. This test is used to determine a prior exposure to VZV and to aid in the determination of acute or convalescent stage of VZV infection. There are several situations in which providing a specific laboratory diagnosis of VZV infection is crucial. The first is in patients who are receiving immunosuppressive therapy or who have abnormalities in their cell-mediated immune responses and the second is in children receiving chemotherapy and radiotherapy for cancer. VZV infection in these cases may cause severe disease or be fatal. Providing a specific diagnosis of VZV infection in these cases may guide in the administration of anti-viral agents or other treatments.

SPECIMEN REQUIREMENTS:

2ml collected in a serum separator tube (gel barrier). Separate serum from cells ASAP or within 2 hours of collection by centrifugation. Stability after separation from cells: Ambient: 48 Hours; Refrigerated: 1 week; Frozen: 1 year (avoid repeated freeze/thaw cycles).

REJECTION CRITERIA: Plasma or other body fluids. Gross hemolysis

METHOD: ELISA

REFERENCES:

1. Gershon, A.A., LaRussa, P. and Steinberg, S. P. 1995. Varicella-Zoster Virus. In: Manual of Clinical Microbiology. Murray, P.R., Baron, E.J., Tenover, F.C. and Tenover, R.H. (eds). 6th Edition, ASM Press, Washington, DC. p.895-904.
2. Heath R. B. 1987. Varicella-Zoster. In: Principle and Practice of Clinical Virology. Zuckerman, A.J., Banatvala, J. E. and Pattison, J.R. (eds) John Wiley and Sons Ltd., New York, p51-73.
3. Arvin, A. M., Koropchak, C.M. and Wittek, A. E. 1983. Immunologic Evidence of Reinfection with Varicella-Zoster Virus. J. Infect. Dis. 148, No. 2: 200-205.
4. Bio-Rad Measles IgG EIA, www.bio-rad.com/webroot/web/pdf/cdg/literature/J-114A_VZV.pdf

RESULTS AND INTERPRETTION:

Index Value	Interpretation
< 0.9	NEGATIVE for VZV IgG, presumed NON-IMMUNE to VZV infection
≥ 0.9 and < 1.1	EQUIVOCAL. Another specimen should be tested 10 to 14 days later in parallel with the initial specimen. If the second specimen is Equivocal, the individual is negative for primary or recent VZV infection and Equivocal for antibody status. If the 2nd sample is positive, the individual can be considered to have a primary infection.
≥ 1.1	POSITIVE for VZV IgG, presumed IMMUNE to VZV infection

Turnaround time: 7 business days